

# FORENSIC BIOLOGY PROTOCOLS FOR FORENSIC STR TYPING

<b>GENEMAPPER ID – RERUN CODES</b>		
DATE EFFECTIVE 07/11/2016	APPROVED BY NUCLEAR DNA TECHNICAL LEADER	PAGE 1 OF 1

<b>ReRun Codes</b>	
Sample Status	Code
All peak labels removed.	**
Peak(s) within basepair range affected by overblown peak(s). Labels removed.	*
Sample shows presence of OL allele	^
No or poor size standard	#

System for Rerun	Code
Identifiler	I
MiniFiler	F
Yfiler	M
Do not rerun	N/A

Parameter for Rerun	Code
Normal (HCN)	no code
High (HCN)	R
1/5 dilution	D.2
1/10 dilution	D.1
1/20 dilution	D.05
1/100 dilution	D.01
Re-aliquot 1 ul	1ul
Re-aliquot 2 ul	2ul
1 kV 22 s (LCN)	L
3 kV 20 s (LCN)	N
6 kV 30 s (LCN)	H

Controlled versions of Department of Forensic Biology Manuals only exist in the Forensic Biology Qualtrax software. All printed versions are non-controlled copies.